## en English

**PRE-FILLED HUMIDIFIER WITH DISPOSABLE HUMIDIFIER ADAPTER**  
**STERILE WATER FOR INHALATION, USP**

**These instructions contain important information for safe use of the product. Read the entire contents of these Instructions For Use, including Warnings and Cautions, before using this product. Failure to properly follow warnings, cautions and instructions could result in death or serious injury to the patient and/or clinician.**

**NOTE: DISTRIBUTE THIS INSERT SHEET TO ALL PRODUCT LOCATIONS.**

**1. DESCRIPTION:**

Disposable, single patient use Respiratory Therapy Bottle with Humidifier Adapter. The disposable humidifier adapter is provided non-sterile. The internal contents of the bottle are sterile unless opened or damaged.

**2. INDICATIONS FOR USE:**

The product is intended to humidify a patient's oxygen flow during oxygen administration and is intended for single patient use in Respiratory Therapy.

**3. CONTRAINDICATIONS:**

None known.

**4. WARNINGS:**

4.1 These products are intended for inhalation only, not for parenteral use or for injections or in preparations to be used for injection. Incorrect use can result in patient injury.

4.2 This product is intended for single patient use. Discard unused water. Failure to do so may result in patient injury.

4.3 Ensure all connections are secure. Failure to do so may result in inadequate oxygen flow to the patient.

**5. CAUTIONS:**

5.1 Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

5.2 Do not use if the product or container is opened or damaged.

5.3 Avoid excessive heat or freezing as this may damage the bottle and result in leaking.

5.4 Follow standard infection control procedures as specified by the Centers for Disease Control and Prevention (USA) or local equivalent.

5.5 Do Not Reuse: Medical devices require specific material characteristics to perform as intended. These characteristics have been verified for single use only and are not to be re-processed for the device for subsequent re-use may adversely affect the integrity of the device or lead to deterioration in performance.

**6. SETUP AND USE:**

6.1 Snap off the outlet port trigger from the bottle and discard (do not twist off).

6.2 Attach the adapter to the inlet of the bottle. Check for secure fit before attaching the oxygen nut to flow meter.

6.3 Attach to flow meter. Set flow rate to 4 lpm and verify airflow through outlet port of the bottle.

6.4 Confirm the alarm function prior to use. Leaving the flow rate at 4 lpm, occlude the outlet port of the bottle and listen for an audible alarm at sound.

6.5 Attach oxygen tubing to small bore outlet port on bottle.

6.6 Adjust the oxygen flow rate selected by physician.

6.7 Record unit dose and date upon completion of therapy.

**7. STORAGE AND HANDLING:**

7.1 Dispose of this product in a safe manner according to local guidelines for disposal of contaminated waste.

7.2 Store at room temperature. Avoid freezing or excessive heat.

The symbol ® indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries.

**Caution - Do Not Reuse - Not made with natural rubber latex - Do not use if package is damaged - Sterilized using aseptic processing techniques - Do not re-sterilize - Non-pyrogenic - Caution: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician - Catalogue Number - Batch Code - Use by - Quantity - Keep away from sunlight - Keep dry - Manufacturer - Authorized Representative in the European Community**

The symbol ® indicates that this mark is registered with the U.S. Patent and Trademark Office and certain other countries.

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