
AVANOS*

ENTERAL SYRINGE

WITH ENFit® CONNECTOR
FOR ENTERAL USE

Instructions for Use

Rx Only: Federal Law (USA) restricts this device to sale by or on the order of a physician.

Indications for Use

AVANOS® Enteral Syringe with ENFit® Connector is indicated for use as a dispenser, a measuring device, and/or a fluid transfer device. It is used to enterally deliver/remove fluids into or from the body. It is intended to be used in clinical or home care settings by users ranging from clinicians to laypersons (under the supervision of a clinician) in all age groups.

⚠ Warnings

- To ensure dose accuracy, for syringe sizes 5mL and below, the moat of the syringe tip must be free of fluids prior to patient administration (Fig. 1).
- Do not autoclave or re-use AVANOS® Enteral Syringe as this may increase the risk of leakage or biological infection.
- Syringe caps (if used) are a choking hazard. Keep out of reach of children.
- Do not modify syringes as this may lead to leakage, inadequate nutrition delivery or patient harm.
- Do not use syringe if syringe tip is visibly blocked as this may lead to inadequate nutrition or medication delivery.
- Do not use if syringe is damaged as this may lead to fluid contamination or inadequate nutrition/medication delivery.
- Sterile syringes are sterilized using ethylene oxide. Do not resterilize.
- To minimize potential for clogging of tip, mix any supplements or additives thoroughly to achieve uniform consistency.
- Be aware that using a medicine cup to fill may cause fluid or medications to enter the moat of the syringe and lead to possible overdose.
- Use a new syringe to flush the medication or fluid after administering any medication to prevent overdose due to the dead space (the remaining fluid in the syringe after administration) in the syringe.

Contraindications

The AVANOS® Enteral Syringe is intended to only connect to products with ENFit® male connectors.

Do not use this device with small bore connectors from other healthcare applications as it may have the potential to misconnect.

Instructions

1. Follow physician's instructions or hospital protocol in opening package and use of contents.
2. Always use technique that minimizes risk of infection. Syringe should be connected to a male ENFit® connector during filling and removal of air bubbles to reduce the risk of dosing volume inaccuracies.
3. Fully depress the plunger before filling syringe. Draw fluid (e.g., formula or medication) into syringe. For accurate measurement, use the black leading edge of the gasket toward the plunger tip. Do not fill and attempt delivery past the maximum graduation mark.

4. To expel excess air bubbles, hold syringe with tip pointing up. Tap side of syringe to allow bubbles to rise toward tip. Draw back slightly on plunger, then push plunger upward to eject air in barrel and tip. Repeat until bubbles are eliminated and verify dose once more.
5. For syringe sizes 5mL and below, ensure the moat of the syringe tip is free of fluids (Fig. 1) and verify dose (repeat as necessary).

Note: The use of bottle adaptors or straws will improve dose accuracy and air removal.


Fig. 1

6. Low Dose Tip syringes should be tapped, or flicked to remove fluid that may be outside the fluid pathway.
7. If needed, secure the syringe cap (sold separately) on the syringe and label until ready to use (label not included).
8. Before delivery, to minimize potential for clogging, mix any supplements or additives thoroughly to achieve uniform consistency.
9. Dispose of syringe as required by hospital protocols and local or national regulations.

Sterile		Bulk Non-Sterile	
Size	Code Number	Size	Code Number
1 mL	SYR-01S	1 mL	SYR-01NS
3 mL	SYR-03S	3 mL	SYR-03NS
5 mL	SYR-05S	5 mL	SYR-05NS
10 mL	SYR-10S	10 mL	SYR-10NS
20 mL	SYR-20S	20 mL	SYR-20NS
60 mL	SYR-60S	60 mL	SYR-60NS

Single Use Only	Non sterile	Sterilized using ethylene oxide
Do not use if package is damaged	Do not resterilize	Not made with natural rubber latex
Product is not made with DEHP as a plasticizer	Product is not made with BPA	Rx Only
Caution	Consult instructions for use	Keep dry
Keep away from sunlight	Manufacturer	Authorized representative in the european community
Catalogue number	Unique Device Identifier	Batch code
		CE (Conformité Européenne) EC Directive Compliance
		Use-by date

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🏢 Jiangsu Caina Medical Co., Ltd., No.23, Huanxi Rd., Zhutang Town,
Distributed in the USA by Avanos Medical Sales, LLC, 5405 Windward Parkway,
[EC REP] Shanghai International Holding Corp. GmbH (Europe),