




Suction Catheters

Mully with Funnel



Area of application

Suction catheters are used for oro-, nasopharyngeal and tracheobronchial suctioning. Must be attached to a jar/canister via a tube before being connected to vacuum source. The products are used for single use only.

Product range

Description	CH	Connector Color	Tip	Effective Length mm	Packaging	Standard (SH A 76) Reference	SOFT (SH A 64) Reference
Suction Catheter Mully with Funnel, Straight tip, Effective length 490mm							
Mully with funnel, straight tip, CH6	6	Light green	Straight	490	1/100/600	50000218	
Mully with funnel, straight tip, CH8	8	Light blue	Straight	490	1/100/600	50000318	
Mully with funnel, straight tip, CH10	10	Black	Straight	490	1/100/600	50000418	50100118
Mully with funnel, straight tip, CH12	12	White	Straight	490	1/100/600	50000518	50100218
Mully with funnel, straight tip, CH14	14	Green	Straight	490	1/100/600	50000618	50100318
Mully with funnel, straight tip, CH16	16	Orange	Straight	490	1/100/600	50000718	50100418
Mully with funnel, straight tip, CH18	18	Red	Straight	490	1/100/600	50000818	
Suction Catheter Mully with Funnel, Curved tip, Effective length 490mm							
Mully with funnel, curved tip, CH8	8	Light blue	Curved	500	1/100/600	80000118	
Mully with funnel, curved tip, CH10	10	Black	Curved	500	1/100/600	80000218	
Mully with funnel, curved tip, CH12	12	White	Curved	500	1/100/600	80000318	
Mully with funnel, curved tip, CH14	14	Green	Curved	500	1/100/600	80000418	
Mully with funnel, curved tip, CH16	16	Orange	Curved	500	1/100/600	80000518	
Mully with funnel, curved tip, CH18	18	Red	Curved	500	1/100/600	80000618	
Suction Catheter Mully with Funnel, Straight tip, Effective length 570mm							
Mully with funnel, straight tip, CH6	6	Light green	Straight	570	1/100/600	50001028	
Mully with funnel, straight tip, CH8	8	Light blue	Straight	570	1/100/600	50001128	
Mully with funnel, straight tip, CH10	10	Black	Straight	570	1/100/600	50001228	50100628
Mully with funnel, straight tip, CH12	12	White	Straight	570	1/100/600	50001328	50100728
Mully with funnel, straight tip, CH14	14	Green	Straight	570	1/100/600	50001428	50100828
Mully with funnel, straight tip, CH16	16	Orange	Straight	570	1/100/600	50001528	50100928
Mully with funnel, straight tip, CH18	18	Red	Straight	570	1/100/600	50001628	50101028
Suction Catheter Mully with Funnel, Curved tip, Effective length 570mm							
Mully with funnel, curved tip, CH10	10	Black	Curved	570	1/100/600	80000728	80200728
Mully with funnel, curved tip, CH12	12	White	Curved	570	1/100/600	80000828	80200828
Mully with funnel, curved tip, CH14	14	Green	Curved	570	1/100/600	80000928	80200928
Mully with funnel, curved tip, CH16	16	Orange	Curved	570	1/100/600	80001028	80201028
Mully with funnel, curved tip, CH18	18	Red	Curved	570	1/100/600	80001128	



Suction Catheters

Mully with Funnel

Product range continued

Description	CH	Connector Color	Tip	Effective Length mm	Packaging	Standard (SH A 76) Reference	SOFT (SH A 64) Reference
Suction Catheter Mully Metric with Funnel, Straight tip, Effective length 330mm, with 1-cm (METRIC) markings							
Mully Metric with funnel, straight tip, CH6	6	Light green	Straight	330	1/100/600	12000328	
Mully Metric with funnel, straight tip, CH8	8	Light blue	Straight	330	1/100/600	12000428	
Mully Metric with funnel, straight tip, CH10	10	Black	Straight	330	1/100/600	12000528	12400128
Mully Metric with funnel, straight tip, CH12	12	White	Straight	330	1/100/600		12400228
Mully Metric with funnel, straight tip, CH14	14	Green	Straight	330	1/100/600		12400328
Mully Metric with funnel, straight tip, CH16	16	Orange	Straight	330	1/100/600		12400428
Suction Catheter Mully with Funnel, Straight tip, Effective length 210mm (Tracheostomy length)							
Mully with funnel, straight tip, CH12	12	White	Straight	210	1/100/600	50020128	
Mully with funnel, straight tip, CH14	14	Green	Straight	210	1/100/600	50020228	

Sterilization

Sterile, Sterilized by ethylene oxide, Shelf life: 5 years

Certifications

CE-marking	CE 0123
Product classification	Class IIA (EU)
Quality system certified to	ISO 13485, Medical Device Single Audit Program (MDSAP)
European Medical Device Nomenclature (EMDN)	A06010103 CANNULAS AND TIPS FOR APIRATION AND IRRIGATION, SINGLE-USE
Global Medical Device Nomenclature (GMDN)	10749

Instructions for use

Included?	Yes
-----------	-----

Material specifications

Contains latex?	No
Contains phthalates?	No
Contains bisphenols?	No
Contains cotton?	No
Contains paper?	Packaging only
Tube and Connector	Polyvinyl chloride

Packaging material

Primary package	Coated paper, polypropylene/polyethylene film
Secondary package	Single-wall corrugated fiberboard
Tertiary package	Double-wall corrugated fiberboard

Suction Catheters

Mully with Funnel



MEDICAL EQUIPMENT

Packaging information									
REF	Primary pack./Single peelpack		Secondary pack./Shelf box		Tertiary pack./Shipper				Pallet Qty. pcs
	Qty. pcs	Dimensions, mm (LxW)	Qty. pcs	Dimensions, mm, (LxWxD)	Qty. pcs	Dimensions mm (LxWxD)	Volume m³	Gross weight kg	
Suction Catheter Mully with Funnel, Straight tip, Effective length 490mm									
50000218	1	630x40	100	640x175x120	600	660x375x370	0.092	7.0	9,000
50000318	1	630x40	100	640x175x120	600	660x375x370	0.092	7.5	9,000
50000418	1	630x40	100	640x175x120	600	660x375x370	0.092	8.0	9,000
50000518	1	630x40	100	640x175x120	600	660x375x370	0.092	8.5	9,000
50000618	1	630x40	100	640x175x120	600	660x375x370	0.092	9.5	9,000
50000718	1	630x40	100	640x175x120	600	660x375x370	0.092	10.5	9,000
50000818	1	630x40	100	640x175x120	600	660x375x370	0.092	11.5	9,000
50100118	1	630x40	100	640x175x120	600	660x375x370	0.092	8.0	9,000
50100218	1	630x40	100	640x175x120	600	660x375x370	0.092	8.5	9,000
50100318	1	630x40	100	640x175x120	600	660x375x370	0.092	9.5	9,000
50100418	1	630x40	100	640x175x120	600	660x375x370	0.092	10.5	9,000
Suction Catheter Mully with Funnel, Curved tip, Effective length 490mm									
80000118	1	630x40	100	640x175x120	600	660x375x370	0.092	7.5	9,000
80000218	1	630x40	100	640x175x120	600	660x375x370	0.092	8.0	9,000
80000318	1	630x40	100	640x175x120	600	660x375x370	0.092	8.5	9,000
80000418	1	630x40	100	640x175x120	600	660x375x370	0.092	9.5	9,000
80000518	1	630x40	100	640x175x120	600	660x375x370	0.092	10.5	9,000
80000618	1	630x40	100	640x175x120	600	660x375x370	0.092	11.5	9,000
Suction Catheter Mully with Funnel, Straight tip, Effective length 570mm									
50001028	1	713x40	100	720x175x120	600	740x375x370	0.103	7.5	9,000
50001128	1	713x40	100	720x175x120	600	740x375x370	0.103	7.9	9,000
50001228	1	713x40	100	720x175x120	600	740x375x370	0.103	8.6	9,000
50001328	1	713x40	100	720x175x120	600	740x375x370	0.103	9.6	9,000
50001428	1	713x40	100	720x175x120	600	740x375x370	0.103	10.4	9,000
50001528	1	713x40	100	720x175x120	600	740x375x370	0.103	11.8	9,000
50001628	1	713x40	100	720x175x120	600	740x375x370	0.103	12.8	9,000
50100628	1	713x40	100	720x175x120	600	740x375x370	0.103	8.6	9,000
50100728	1	713x40	100	720x175x120	600	740x375x370	0.103	9.6	9,000
50100828	1	713x40	100	720x175x120	600	740x375x370	0.103	10.4	9,000
50100928	1	713x40	100	720x175x120	600	740x375x370	0.103	11.8	9,000
50101028	1	713x40	100	720x175x120	600	740x375x370	0.103	12.8	9,000
Suction Catheter Mully with Funnel, Curved tip, Effective length 570mm									
80000728	1	713x40	100	720x175x120	600	740x375x370	0.103	8.6	9,000
80000828	1	713x40	100	720x175x120	600	740x375x370	0.103	9.6	9,000
80000928	1	713x40	100	720x175x120	600	740x375x370	0.103	10.4	9,000
80001028	1	713x40	100	720x175x120	600	740x375x370	0.103	11.8	9,000
80001128	1	713x40	100	720x175x120	600	740x375x370	0.103	12.8	9,000
80200728	1	713x40	100	720x175x120	600	740x375x370	0.103	8.6	9,000
80200828	1	713x40	100	720x175x120	600	740x375x370	0.103	9.6	9,000
80200928	1	713x40	100	720x175x120	600	740x375x370	0.103	10.4	9,000
80201028	1	713x40	100	720x175x120	600	740x375x370	0.103	11.8	9,000

Suction Catheters

Mully with Funnel

SINGLE - USE



Packaging information continued

REF	Primary pack./Single peelpack		Secondary pack./Shelf box		Tertiary pack./Shipper				Pallet Qty. pcs
	Qty. pcs	Dimensions, mm (LxW)	Qty. pcs	Dimensions, mm, (LxWxD)	Qty. pcs	Dimensions mm (LxWxD)	Volume m³	Gross weight kg	
Suction Catheter Mully Metric with Funnel, Straight tip, Effective length 330mm, with 1-cm (METRIC) markings									
12000328	1	485×33.5	100	500×175×120	600	520×375×370	0.072	6.5	12,000
12000428	1	485×33.5	100	500×175×120	600	520×375×370	0.072	7.0	12,000
12000528	1	485×33.5	100	500×175×120	600	520×375×370	0.072	7.5	12,000
12400128	1	485×33.5	100	500×175×120	600	520×375×370	0.072	7.5	12,000
12400228	1	485×33.5	100	500×175×120	600	520×375×370	0.072	8.3	12,000
12400328	1	485×33.5	100	500×175×120	600	520×375×370	0.072	8.7	12,000
12400428	1	485×33.5	100	500×175×120	600	520×375×370	0.072	9.5	12,000
Suction Catheter Mully with Funnel, Straight tip, Effective length 210mm (Tracheostomy length)									
50020128	1	322×33.5	100	325×175×120	600	380×375×370	0.053	7.5	18,000
50020228	1	322×33.5	100	325×175×120	600	380×375×370	0.053	8.0	18,000

GTIN Codes

REF	Primary package	Secondary package	Tertiary package
Suction Catheter Mully with Funnel, Straight tip, Effective length 490mm			
50000218	06944932786523	16944932786520	26944932786527
50000318	06944932786530	16944932786537	26944932786534
50000418	06944932786547	16944932786544	26944932786541
50000518	06944932786554	16944932786551	26944932786558
50000618	06944932786561	16944932786568	26944932786565
50000718	06944932786578	16944932786575	26944932786572
50000818	06944932786585	16944932786582	26944932786589
50100118	06936881153990	16936881153997	26936881153994
50100218	06936881154003	16936881154000	26936881154007
50100318	06936881154010	16936881154017	26936881154014
50100418	06936881154027	16936881154024	26936881154021
Suction Catheter Mully with Funnel, Curved tip, Effective length 490mm			
80000118	06936881154195	16936881154192	26936881154199
80000218	06936881154201	16936881154208	26936881154205
80000318	06936881154218	16936881154215	26936881154212
80000418	06936881154225	16936881154222	26936881154229
80000518	06936881154232	16936881154239	26936881154236
80000618	06936881154249	16936881154246	26936881154243

Suction Catheters

Mully with Funnel



GTIN Codes continued			
REF	Primary package	Secondary package	Tertiary package
Suction Catheter Mully with Funnel, Straight tip, Effective length 570mm			
50001028	6936881159473	16936881159470	26936881159477
50001128	6936881159480	16936881159487	26936881159484
50001228	6936881159497	16936881159494	26936881159491
50001328	6936881159503	16936881159500	26936881159507
50001428	6936881159510	16936881159517	26936881159514
50001528	6936881159527	16936881159524	26936881159521
50001628	6936881159534	16936881159531	26936881159538
50100628	6936881159541	16936881159548	26936881159545
50100728	6936881159558	16936881159555	26936881159552
50100828	6936881159565	16936881159562	26936881159569
50100928	6936881159572	16936881159579	26936881159576
50101028	6936881159589	16936881159586	26936881159583
Suction Catheter Mully with Funnel, Curved tip, Effective length 570mm			
80000728	6936881159749	16936881159746	26936881159743
80000828	6936881159756	16936881159753	26936881159750
80000928	6936881159763	16936881159760	26936881159767
80001028	6936881159770	16936881159777	26936881159774
80001128	6936881159787	16936881159784	26936881159781
80200728	6936881159701	16936881159708	26936881159705
80200828	6936881159718	16936881159715	26936881159712
80200928	6936881159725	16936881159722	26936881159729
80201028	6936881159732	16936881159739	26936881159736
Suction Catheter Mully with Funnel, Straight tip, Effective length 330mm, with 1-cm (METRIC) markings			
12000328	6936881160011	16936881160018	26936881160015
12000428	6936881160028	16936881160025	26936881160022
12000528	6936881160035	16936881160032	26936881160039
12400128	6936881160042	16936881160049	26936881160046
12400228	6936881160059	16936881160056	26936881160053
12400328	6936881160066	16936881160063	26936881160060
12400428	6936881160073	16936881160070	26936881160077
Suction Catheter Mully with Funnel, Straight tip, Effective length 210mm (Tracheostomy length)			
50020128	6936881160165	16936881160162	26936881160169
50020228	6936881160172	16936881160179	26936881160176

Suction Catheters

Mully with Funnel

SINGLE - USE



MEDICAL EQUIPMENT

Storage and handling	
	THIS WAY UP
	Keep away from rain
	Fragile; handle with care
	Keep away from sunlight
	Do not resterilize
	Do not re-use
	Do not use if package is damaged and consult electronic instructions for use
	Temperature limit (-15° to 49°C)
	Consult instructions for use or use electronic instructions for use

SPONSORED IN AUSTRALIA BY

Unoquip Pty Ltd
PO Box 1117
Windsor VIC 3181

1300 55 80 79
info@unoquip.au

DISTRIBUTED BY

UNOQUIP GMBH
Schützengraben 20
8200 Schaffhausen
Switzerland
info@unoquip.com
www.unoquip.com



Well Lead Medical Co., Ltd.
C-4 Jinhu Industrial Estate
Hualong, Panyu Guangzhou
511434 P.R. China



Shanghai Intl. Holding
Corp. GmbH
Eiffestrasse 80
20537 Hamburg
Germany



QS Engineering AG
Erlenstrasse 31
CH - 4106 Therwil
Switzerland

**UK
Responsible
Person
(UKRP)**

SUNGO Certification
Company Limited
3rd Floor, 70 Gracechurch
Street, London, EC3V 274HR
England, United Kingdom