


Intermittent Catheters



Area of application

Intermittent Catheters are used to provide an intermittent pathway for short term drainage of urine from the bladder of male (Nelaton and Tiemann) and female (Female) patients. The products are for single use only.

Product range

Size	Nelaton Catheters		Nelaton Catheters with Sheath		Tiemann Catheters	Female catheters
	Standard (Shore A 76)	SOFT (Shore A 64)	Standard (Shore A 76)	SOFT (Shore A 64)	Standard (Shore A 76)	Standard (Shore A 76)
CH/FG	Length 400mm	Length 400mm	Length 400mm	Length 400mm	Length 400mm	Length 225mm
6	11300118		11400118		11700118	11200118
8	11300218		11400218		11700218	11200218
10	11300318	11600318	11400318	11500318	11700318	11200318
12	11300418	11600418	11400418	11500418	11700418	11200418
14	11300518	11600518	11400518	11500518	11700518	11200518
16	11300618	11600618	11400618	11500618	11700618	11200618
18	11300718	11600718	11400718	11500718	11700718	11200718
20	11300818		11400818		11700818	

Sterilization

Sterile, Sterilized by ethylene oxide, Shelf life: 5 years

Certifications

CE-marking	CE 0123
Product classification	Class IIA (EU)
Quality system certified to	ISO 13485 Medical Device Single Audit Program (MDSAP)
European Medical Device Nomenclature (EMDN)	U01010502 - NELATON CATHETERS, NON SELF-LUBRICATING
Global Medical Device Nomenclature (GMDN)	34930

Instructions for use

Included?	Yes
-----------	-----

Material specifications

Contains latex?	No
Contains DEHP?	No
Contains cotton?	No
Contains paper?	Packaging only
Catheter shaft and funnel	Polyvinyl chloride
Sheath	Polyethylene



Intermittent Catheters

Packaging information

Reference	Primary pack./Single peelpack		Secondary pack./Shelf box		Tertiary pack./Shipper			Gross weight tkg	Pallet Qty. pcs
	Qty. pcs	Dimensions, mm (LxW)	Qty. pcs	Dimensions mm, (LxWxD)	Qty. pcs	Dimensions mm (LxWxD)	Volume m³		
11300118	1	485 × 33.5	100	500 × 175 × 120	600	520 × 375 × 370	0.0722	5.0	12,000
11300218	1	485 × 33.5	100	500 × 175 × 120	600	520 × 375 × 370	0.0722	5.5	12,000
11300318	1	485 × 33.5	100	500 × 175 × 120	600	520 × 375 × 370	0.0722	6.0	12,000
11300418	1	485 × 33.5	100	500 × 175 × 120	600	520 × 375 × 370	0.0722	6.5	12,000
11300518	1	485 × 33.5	100	500 × 175 × 120	600	520 × 375 × 370	0.0722	7.0	12,000
11300618	1	485 × 33.5	100	500 × 175 × 120	600	520 × 375 × 370	0.0722	8.0	12,000
11300718	1	485 × 33.5	100	500 × 175 × 120	600	520 × 375 × 370	0.0722	9.0	12,000
11300818	1	485 × 33.5	100	500 × 175 × 120	600	520 × 375 × 370	0.0722	9.5	12,000
11600318	1	485 × 33.5	100	500 × 175 × 120	600	520 × 375 × 370	0.0722	9.5	12,000
11600418	1	485 × 33.5	100	500 × 175 × 120	600	520 × 375 × 370	0.0722	9.5	12,000
11600518	1	485 × 33.5	100	500 × 175 × 120	600	520 × 375 × 370	0.0722	9.5	12,000
11600618	1	485 × 33.5	100	500 × 175 × 120	600	520 × 375 × 370	0.0722	9.5	12,000
11600718	1	485 × 33.5	100	500 × 175 × 120	600	520 × 375 × 370	0.0722	9.5	12,000
11400118	1	485 × 33.5	100	500 × 175 × 120	600	520 × 375 × 370	0.0722	9.5	12,000
11400218	1	485 × 33.5	100	500 × 175 × 120	600	520 × 375 × 370	0.0722	9.5	12,000
11400318	1	485 × 33.5	100	500 × 175 × 120	600	520 × 375 × 370	0.0722	9.5	12,000
11400418	1	485 × 33.5	100	500 × 175 × 120	600	520 × 375 × 370	0.0722	9.5	12,000
11400518	1	485 × 33.5	100	500 × 175 × 120	600	520 × 375 × 370	0.0722	9.5	12,000
11400618	1	485 × 33.5	100	500 × 175 × 120	600	520 × 375 × 370	0.0722	9.5	12,000
11400718	1	485 × 33.5	100	500 × 175 × 120	600	520 × 375 × 370	0.0722	9.5	12,000
11400818	1	485 × 33.5	100	500 × 175 × 120	600	520 × 375 × 370	0.0722	9.5	12,000
11500318	1	485 × 33.5	100	500 × 175 × 120	600	520 × 375 × 370	0.0722	9.5	12,000
11500418	1	485 × 33.5	100	500 × 175 × 120	600	520 × 375 × 370	0.0722	9.5	12,000
11500518	1	485 × 33.5	100	500 × 175 × 120	600	520 × 375 × 370	0.0722	9.5	12,000
11500618	1	485 × 33.5	100	500 × 175 × 120	600	520 × 375 × 370	0.0722	9.5	12,000
11500718	1	485 × 33.5	100	500 × 175 × 120	600	520 × 375 × 370	0.0722	9.5	12,000
11700118	1	485 × 33.5	100	500 × 175 × 120	600	520 × 375 × 370	0.0722	4.0	12,000
11700218	1	485 × 33.5	100	500 × 175 × 120	600	520 × 375 × 370	0.0722	4.0	12,000
11700318	1	485 × 33.5	100	500 × 175 × 120	600	520 × 375 × 370	0.0722	4.5	12,000
11700418	1	485 × 33.5	100	500 × 175 × 120	600	520 × 375 × 370	0.0722	4.5	12,000
11700518	1	485 × 33.5	100	500 × 175 × 120	600	520 × 375 × 370	0.0722	5.0	12,000
11700618	1	485 × 33.5	100	500 × 175 × 120	600	520 × 375 × 370	0.0722	5.5	12,000
11700718	1	485 × 33.5	100	500 × 175 × 120	600	520 × 375 × 370	0.0722	6.0	12,000
11700818	1	485 × 33.5	100	500 × 175 × 120	600	520 × 375 × 370	0.0722	6.5	12,000
11200118	1	322 × 33.5	100	325 × 175 × 120	600	380 × 345 × 370	0.0485	4.0	18,000
11200218	1	322 × 33.5	100	325 × 175 × 120	600	380 × 345 × 370	0.0485	4.0	18,000
11200318	1	322 × 33.5	100	325 × 175 × 120	600	380 × 345 × 370	0.0485	4.5	18,000
11200418	1	322 × 33.5	100	325 × 175 × 120	600	380 × 345 × 370	0.0485	4.5	18,000
11200518	1	322 × 33.5	100	325 × 175 × 120	600	380 × 345 × 370	0.0485	5.0	18,000
11200618	1	322 × 33.5	100	325 × 175 × 120	600	380 × 345 × 370	0.0485	5.5	18,000
11200718	1	322 × 33.5	100	325 × 175 × 120	600	380 × 345 × 370	0.0485	6.0	18,000

Intermittent Catheters

GTIN codes			
REF	Primary package	Secondary package	Tertiary package
11300118	06936881134517	16936881134514	26936881134511
11300218	06936881134524	16936881134521	26936881134528
11300318	06936881134531	16936881134538	26936881134535
11300418	06936881134548	16936881134545	26936881134542
11300518	06936881134555	16936881134552	26936881134559
11300618	06936881134562	16936881134569	26936881134566
11300718	06936881134579	16936881134576	26936881134573
11300818	06936881134586	16936881134583	26936881134580
11600318	06936881152979	16936881152976	26936881152973
11600418	06936881152986	16936881152983	26936881152980
11600518	06936881152993	16936881152990	26936881152997
11600618	06936881153006	16936881153003	26936881153000
11600718	06936881153013	16936881153010	26936881153017
11400118	06936881152849	16936881152846	26936881152843
11400218	06936881152856	16936881152853	26936881152850
11400318	06936881152863	16936881152860	26936881152867
11400418	06936881152870	16936881152877	26936881152874
11400518	06936881152887	16936881152884	26936881152881
11400618	06936881152894	16936881152891	26936881152898
11400718	06936881152900	16936881152907	26936881152904
11400818	06936881152917	16936881152914	26936881152911
11500318	06936881152924	16936881152921	26936881152928
11500418	06936881152931	16936881152938	26936881152935
11500518	06936881152948	16936881152945	26936881152942
11500618	06936881152955	16936881152952	26936881152959
11500718	06936881152962	16936881152969	26936881152966
11700118	06936881134616	16936881134613	26936881134610
11700218	06936881134623	16936881134620	26936881134627
11700318	06936881134630	16936881134637	26936881134634
11700418	06936881134647	16936881134644	26936881134641
11700518	06936881134654	16936881134651	26936881134658
11700618	06936881134661	16936881134668	26936881134665
11700718	06936881134678	16936881134675	26936881134672
11700818	06936881134685	16936881134682	26936881134689
11200118	06936881134418	16936881134415	26936881134412
11200218	06936881134425	16936881134422	26936881134429
11200318	06936881134432	16936881134439	26936881134436
11200418	06936881134449	16936881134446	26936881134443
11200518	06936881134456	16936881134453	26936881134450
11200618	06936881134463	16936881134460	26936881134467
11200718	06936881134470	16936881134477	26936881134474

Packaging material	
Primary package	Coated paper, Film
Secondary package	Ivory Cardboard
Tertiary package	Double-wall corrugated fiberboard

Intermittent Catheters

Storage and handling	
	This way up
	Keep dry
	Fragile, handle with care
	Keep away from sunlight
	Do not resterilize
	Do not re-use
	Do not use if package is damaged and consult instructions for use
	Temperature limit (-15° to 49°C)
	Consult instructions for use or consult electronic instructions for use

SPONSORED IN AUSTRALIA BY

Unoquip Pty Ltd
PO Box 1117
Windsor VIC 3181


1300 55 80 79
info@unoquip.au

DISTRIBUTED BY

UNOQUIP GMBH
Schützengraben 20
8200 Schaffhausen
Switzerland
info@unoquip.com
www.unoquip.com



Well Lead Medical Co., Ltd.
C-4 Jinhu Industrial Estate
Hualong, Panyu Guangzhou
511434 P.R. China



Shanghai Intl. Holding
Corp. GmbH
Eiffestrasse 80
20537 Hamburg
Germany



QS Engineering AG
Erlenstrasse 31
CH - 4106 Therwil
Switzerland

**UK
Responsible
Person
(UKRP)**

SUNGO Certification
Company Limited
3rd Floor, 70 Gracechurch
Street, London, EC3V 274HR
England, United Kingdom