




Duodenal Tubes



Area of application

Duodenal tubes are used to drain gastric contents, decompress the stomach or duodenum, obtain a specimen of the gastric contents and administer medicine, liquids and/or enteral feeding. They are inserted through the nose or mouth via the oesophagus into the stomach or duodenum. The product is for single use only.

Product range

Tube Outer Diameter, mm	CH/FG Size Color-coded funnel connector	Levin	Ryles
		X-ray line 2 eyes open end Length: 1,200mm*	X-ray line and tip, Luer adapter 4 eyes closed end Length: 1,210mm*
		Reference	Reference
Duodenal Tube with graduation marks, Shore A 84			
2.67	8	10000128	11000128
3.33	10	10000228	11000228
4.00	12	10000328	11000328
4.67	14	10000428	11000428
5.33	16	10000528	11000528
6.00	18	10000628	11000628

*Total length of tube including connector

Sterilization

Sterile, Sterilized by ethylene oxide, Shelf life: 5 years

Certifications

CE-marking	CE 0123
Product classification	Class IIA (EU)
Quality system certified to	ISO 13485
European Medical Device Nomenclature (EMDN)	G020299
Global Medical Device Nomenclature (GMDN)	14221

Instructions for use

Included?	Yes
-----------	-----

Material specifications

Contains latex?	No
Contains DEHP?	No
Contains cotton?	No
Contains paper?	Yes (packaging only)
Tube	Polyvinyl chloride
Connector	Polyvinyl chloride



Duodenal Tubes

Packaging information

Reference	Primary pack./Single peelpack		Secondary pack./Shelf box		Tertiary pack./Shipper			Gross weight kg	Pallet Qty. pcs
	Qty. pcs	Dimensions, mm (LxW)	Qty. pcs	Dimensions mm (LxWxD)	Qty. pcs	Dimensions mm (LxWxD)	Volume m ³		
10000128	1	230x130	50	300x230x130	200	540x250x320	0.043	3.73	7,200
10000228	1	230x130	50	300x230x130	200	540x250x320	0.043	4.19	7,200
10000328	1	230x130	50	300x230x130	200	540x250x350	0.047	4.83	6,000
10000428	1	230x130	30	300x230x130	120	540x250x320	0.043	3.85	4,320
10000528	1	230x130	30	300x230x130	120	540x250x320	0.043	4.23	4,320
10000628	1	230x130	30	300x230x130	120	540x250x320	0.043	4.74	4,320
11000128	1	230x130	50	300x230x130	200	540x250x320	0.043	4.53	7,200
11000228	1	230x130	50	300x230x130	200	540x250x320	0.043	4.99	7,200
11000328	1	230x130	50	300x230x130	200	540x250x350	0.047	5,63	6,000
11000428	1	230x130	30	300x230x130	120	540x250x320	0.043	4,33	4,320
11000528	1	230x130	30	300x230x130	120	540x250x320	0.043	4.71	4,320
11000628	1	230x130	30	300x230x130	120	540x250x320	0.043	5.22	4,320

GTIN codes

REF	Primary package	Secondary package	Tertiary package
10000128	16970754490864	26970754490861	36970754490868
10000228	16970754490871	26970754490878	36970754490875
10000328	16970754490888	26970754490885	36970754490882
10000428	16970754490895	26970754490892	36970754490899
10000528	16970754490901	26970754490908	36970754490905
10000628	16970754490918	26970754490915	36970754490912
11000128	16970754490925	26970754490922	36970754490929
11000228	16970754490932	26970754490939	36970754490936
11000328	16970754490949	26970754490946	36970754490943
11000428	16970754490956	26970754490953	36970754490950
11000528	16970754490963	26970754490960	36970754490967
11000628	16970754490970	26970754490977	36970754490974

Packaging material

Primary package	Arjor paper, Polyethylene terephthalate/Polyethylene film
Secondary package	Ivory cardboard
Tertiary package	Corrugated cardboard

Duodenal Tubes



Storage and handling	
	Keep away from rain
	Keep away from sunlight
	Do not re-sterilize
	Do not re-use
	Do not use if package is damaged and consult instructions for use
	Consult instructions for use or consult electronic instructions for use

SPONSORED IN AUSTRALIA BY

Unoquip Pty Ltd
PO Box 1117
Windsor VIC 3181

1300 55 80 79
info@unoquip.au

DISTRIBUTED BY

UNOQUIP GMBH
Schützengraben 20
8200 Schaffhausen
Switzerland
info@unoquip.com
www.unoquip.com



GCMEDICA ENTERPRISE
LTD.(WUXI)
Loujin Industrial Park,
Shuofang, Wuxi,
214143, Jiangsu, P.R. China



Kingsmead Service B.V.
Zonnehof 36, 2632 BE,
Nootdorp, Netherland



CMC Medical Devices GmbH
Bahnhofstrasse 32
CH-6300 Zug
Switzerland

UK
Responsible
Person
(UKRP)

Kingsmead Service Limited
19 Mezzanine Floor 19-21
Crawford Street London
England W1H 1PJ