


PRODUCT SPECIFICATION		This document is the property of Tytex. The contents must not be disclosed to third parties without the expressed written consent from Tytex. Infringements in the above will result in prosecution.	
Product Name: Care Short		Article no: 3022	Article version: 53
APPROVED BY TYTEX		APPROVED BY CUSTOMER	VALID FROM 12. January 2020
		VERSION NO. 1.9	

1. SCOPE

Product to be worn as part of a two-piece reusable or disposable absorbent incontinence solution.

2. INTENDED USE

Textile fixation devices for fixating of i.e. a bandage or a pad, worn as normal underwear pants for alleviation of continence diseases or handi-cap.



Single Patient - multiple use

3. PICTURE



4. MATERIAL COMPOSITION

95 % Polyamide, 5 % Elastane.

Products do not contain latex or natural rubber or pose any harmful human impact associated with its intended use.


Products do not contain human or animal tissue or extracts, prohibited azo-dye pigments, medical substances or phthalates.









The products are in compliance with REACH.

No chemicals added during treatment/finishing. Products may contain softening and antistatic agents up to 1.5 %.

The products are covered by Oeko-Tex Standard 100, class II



PRODUCT SPECIFICATION		This document is the property of Tytex. The contents must not be disclosed to third parties without the expressed written consent from Tytex. Infringements in the above will result in prosecution.	
Product Name: Care Short		Article no: 3022	Article version: 53
APPROVED BY TYTEX		APPROVED BY CUSTOMER	VALID FROM 12. January 2020
		VERSION NO. 1.9	

5. RECOMMENDED CARE Washing at 60° - 95° C/ 140° - 200° F. in washing bag followed by gentle tumble-drying. Recommended symbols (EN 3758):      Best durability achieved at 60° C / 140° F washing. Use of chlorine based bleach or chlorine based disinfection will decrease durability		
6. EXPIRATION Products have an expiry date of five years from the date of manufacturing		
7. TRANSPORTATION, STORAGE AND WASTE HANDLING Transportation and Storage: No special requirements to temperature or humidity. Preferable transport in sealed boxes on pallets with stretch wrapping. Do not stack pallets on top of each other unless previously agreed. Keep the products in a dry, well ventilated place and not exposed to direct sunlight, especially when original stretch wrapping has been removed from packaging. Waste: Products may be disposed with regular household waste		 
8. CLASSIFICATION Products are Medical Devices Class I, None Sterile Design and manufacturing in accordance with Medical Device Regulation 2017/745 and FDA requirements.		
GMDN Code:		35120
FDA Code and Listing:		EYQ E434956

PRODUCT SPECIFICATION

This document is the property of Tytex. The contents must not be disclosed to third parties without the expressed written consent from Tytex. Infringements in the above will result in prosecution.

Product Name: Care Short

Article no: 3022

Article version: 53

APPROVED BY
TYTEX


APPROVED BY
CUSTOMER

VALID FROM

12. January 2020

VERSION NO.


1.9

9. METHOD OF MEASUREMENT

10. SIZE

Parameter	Size		Size Colour	
Size	00	XS	Red	
	01	S	Yellow	
	03	M	Blue	
	05	L	Brown	
	07	XL	Green	

11. TECHNICAL SPECIFICATIONS


Parameter	Size	Unit	Target	Tolerance	Measurement Method
A Waist, relaxed	00	cm	21,0	± 2,0	Ruler
	01		22,0	± 2,0	
	03		25,0	± 2,0	
	05		28,0	± 2,0	
	07		32,0	± 2,0	
B Height, relaxed	00	cm	22,0	± 2,0	Ruler
	01		24,0	± 2,0	
	03		30,0	± 2,0	
	05		34,0	± 2,0	
	07		36,0	± 2,0	
D Crotch width	00	cm	4,0	± 2,0	Ruler
	01		7,0	± 2,0	
	03		7,0	± 2,0	
	05		7,0	± 2,0	
	07		7,0	± 2,0	
A Waist, 1800 g load	00	cm	43,0	± 3,0	On frame
	01		48,0	± 3,0	
	03		55,0	± 3,0	
	05		59,5	± 3,0	
	07		64,0	± 3,0	

PRODUCT SPECIFICATION			This document is the property of Tytex. The contents must not be disclosed to third parties without the expressed written consent from Tytex. Infringements in the above will result in prosecution.	
Product Name: Care Short			Article no: 3022	Article version: 53
APPROVED BY TYTEX		APPROVED BY CUSTOMER	VALID FROM 12. January 2020	VERSION NO. 1.9


11. TECHNICAL SPECIFICATIONS

Parameter	Size	Unit	Target	Tolerance	Measurement Method
C Leg, 1000 g load	00	cm	25,0	± 2,0	On frame
	01		31,0	± 2,0	
	03		36,5	± 2,0	
	05		41,0	± 2,0	
	07		44,0	± 2,0	
Weight	00	g	27,5	± 3,5	Scale
	01		30,0	± 3,5	
	03		40,5	± 3,5	
	05		47,5	± 3,5	
	07		61,0	± 3,5	

12. DESIGN PARAMETERS Not part of daily sampling routine for inspection and test.

Parameter	Size	Unit	Target	Tolerance	Test method
A Waist, 500 g load	00	cm	28,0	± 2,0	On frame
	01		30,0	± 2,0	
	03		34,0	± 2,0	
	05		37,0	± 2,0	
	07		41,5	± 2,0	
C Leg, 200 g load	00	cm	17,0	± 2,0	On frame
	01		21,0	± 2,0	
	03		24,0	± 2,0	
	05		27,0	± 2,0	
	07		30,0	± 2,0	
Height, waistband	00	cm	2,7	± 0,5	Ruler
	01		3,0	± 0,5	
	03		3,0	± 0,5	
	05		3,0	± 0,5	
	07		3,3	± 0,5	
Leg displacement In crotch area. (Sewing)	All	cm	Max 0,5		
Bursting Strength	All	N	350	± 20	ASTM D 3787
Shrinkage after 1st washing*	All	%	7	± 2	Domestic washing
* Durability	All	Washes	50	± 5	Domestic washing

* Assuming proper washing temperature, chemicals and cycle times. (The test does not include the wear and tear of the product)

PRODUCT SPECIFICATION			This document is the property of Tytex. The contents must not be disclosed to third parties without the expressed written consent from Tytex. Infringements in the above will result in prosecution.	
Product Name: Care Short		Article no: 3022	Article version: 53	
APPROVED BY TYTEX		APPROVED BY CUSTOMER	VALID FROM 12. January 2020	VERSION NO. 1.9

13. SIZE RECOMMENDATIONS

Size recommendations are to be used as a guide only. Other size parameters on product packaging may be used at customer discretion.

Size	Waist Measurement		Leg Circumference		Hip Measurement		Body weight	
	cm	Inches	cm / max	Inches	cm	Inches	Kg	Lbs.
00	55-75		35-50		60-90			
01	65-85		40-60		70-110			
03	75-105		50-70		80-130			
05	85-115		60-80		90-150			
07	95-125		70-90		100-170			

14. Waist Band

Front outside	
Inside	  